



HART SCIENTIFIC 5901B-G

Triple Point of Water Cell, Mini Glass, 8mm I.D.

Product Highlights

- Hart Model # 5901B-G
- Expanded Uncertainty (k=2): < 0.0002 °C
- Reproducibility: 0.00005 °C
- Dimensions: 30 mm OD; 8 mm ID; 180 mm long
- Immersion Depth (water surface to well bottom): 118mm
- Material: Borosilicate Glass
- Water Source: Ocean
- dDVSMOW: ±20%
- d18OVSMOW: ±3%
- Effect of deviation from VSMOW: ± 14 μK