



METTLER TOLEDO Classic Light (AL-S) Analytical Balance EW-11334-22

- **Classic Light analytical balances are ideal for easy weighing tasks**
- Automatic adjustment with external weight
- Piece counting, percent weighing, dynamic weighing, \pm weighing, and calculation with a custom free factor
- Weighing units include g, kg, mg, ct, lb, oz, ozt, GN, dwt, momme, mesghal, tael, tical, tola, and baht

Classic Light balances are ideal for easy weighing tasks. Each unit comes standard with a versatile RS232 interface, and offers a choice of applications. No unnecessary features to push up the cost, just an affordable price and the quality you naturally expect from Mettler-Toledo. Balances are powered by a line cord or 4-AA batteries for greater portability (order separately).

What's included: AC power adapter, below balance weighing hook, analytical balances include three-door glass draft shield.

Specifications

Capacity	101g
Readability	0.1 mg
Repeatability	0.1mg
Brand	Mettler Toledo
Linearity	0.3mg
Output	built-in RS-232 interface
Power	115 VAC, 50/60 Hz
Type	Analytical
Housing dimensions	9-3/8"H x 14-5/8"H x 13-1/8"D
Calibration	External
Stability filters	Vibration, drafts
Pan size (" dia)	3-1/2"